


Uncompromized *HbA1c* results



sebia

The new language of life

HbA1c by capillary electrophoresis: *clear-cut and precise results*

Ergonomics

The laboratory optimization consultancy service provides effective solutions that fit your needs: increased productivity, TAT management, hands-on time improvement, consolidation and optimized workflow.

Flexibility

The HbA1c test can be added to a wide range of diagnostic tests on the same platform. This provides flexibility to meet the needs of both: small and large laboratories.

High resolution

The clear separation between HbA1c and the most common hemoglobin variants ensures results without analytical interference. Accurate HbA2 measurement contributes to better patient management.

Productivity

The simultaneous analysis of multiple samples supports productivity, which can be adapted to any workload.


Accuracy

Scientific publications and results of EQA programs reflect the excellent performances of the Sebia method, also IFCC and NGSP certified.

Facilitated interpretation

Thanks to the clarity of the electrophoretic profiles, the time required for validation is reduced, even in the presence of hemoglobin variants or thalassemias, easily and clearly detected.

Pattern with β -Thalassemia






Multiparametric instruments *with large range of tests**:

- HbA1c (diagnosis & monitoring)
- Serum & Urine Protein Electrophoresis (SPE & UPE)
- Serum & Urine Immunotyping
- Hemoglobin Disorders Screening
- CDT/CDT_{IFCC}

Scalable solution *adapted to each lab's activity*



Flexible automation and consolidation

Laboratory optimization consultancy service

"Best in class" separative method

* Check with your Sebia representative for test availability
Maximum Capacity or MC designation refers to: Tube Loader CAPILLARYS 3...1249 / Rack Conveyor CAPILLARYS 3...1250

Ordering information

Instruments

• CAPILLARYS 3 TERA	1246
• CAPILLARYS 3 OCTA	1245
• TUBE LOADER CAPILLARYS 3	1249
• RACK CONVEYOR CAPILLARYS 3	1250
• MINICAP FLEX-PIERCING	1232

Controls/Calibrators

• Hb A1c CAPILLARY CALIBRATORS (2)	4755
• Hb A1c CAPILLARY CALIBRATORS (5X2)	4756
• MULTI-SYSTEM Hb A1c CAPILLARY CONTROLS (10 X 2)	4767
• MULTI-SYSTEM Hb A1c CAPILLARY CONTROLS (2)	4768

Accessories

• CAPILLARYS 3 & MC SWITCH RACK Hb A1c	1383
---	------

Consumables

• CAPILLARYS 3 CAPICLEAN	2060
• CAPILLARYS 3 WASH SOLUTION	2062
• CAPI 3 REAGENT CUPS (24 x 14)	2582
• CAPI 3 BINS FOR USED REAGENT CUPS (5)	2581
• MINICAP FLEX-PIERCING CAPICLEAN	2251
• «AUTOMATIC LOW VOLUME» LABELS FOR MINICAP FLEX-PIERCING (20)	9208
• «MANUAL LOW VOLUME» LABELS FOR MINICAP FLEX-PIERCING (20)	9209
• A1c/CE - CAPILLARY BLOOD SAMPLE KIT (A)*	9212
• A1c/CE LOW VOLUME SAMPLE COLLECTION TUBE	9216
• TEST TUBES (200)	9214

Reagents

• CAPI 3 Hb A1c	2515
• MINICAP Hb A1c	2215

* With the CAPILLARYS 3 automated instrument



Scan here
for more
information
about HbA1c

For further information, please find your Sebia local representative to
contact on our Website www.sebia.com/sebia/worldwide-presence/.

Sebia | RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008
Evry cedex | France | Tel.: +33 1 69 89 80 80 | E-mail: sebia@sebia.com | www.sebia.com

©2026 **Sebia** - RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008 Evry cedex | France - Ref: HBA1CB2EN - 01/2025
Non-contractual illustration - Sebia reserves the right to modify at any time and without notice the information contained in this document intended for healthcare
professionals. This document contains information intended for wide distribution and may therefore contain product details or information that is not available or valid in your
country. Carefully read the instructions in the reagent package inserts and instrument manuals.