


CAPILLARYS
DBS **3**

*The perfect match for screening
hemoglobin disorders
in newborns*



sebia

The new language of life

High throughput solution for newborn hemoglobin disorders screening



Newborn Hemoglobin disorders screening has never been so easy...

SAMPLING

PUNCHER

ELUTION

ANALYSIS

RESULT

Main benefits

1
Ease of use,
speed and
**reduced
hands-on time**

2
Automatic result
flags let you **focus
on pathological
samples**

3
Clear detection of
Hemoglobin variants
(S, C, D, E...) as well as
Thalassemia (Hb Bart's)

Key features



High throughput instrument

70 samples/hr



Full traceability

From sample punch to result



High autonomy

8 full 96-well microplates could be analyzed in a single run



Automatic Reagent Control

RFID for all reagents and Hb AF control



Excellent Resolution

High voltage & temperature-controlled migration




Automatic presumptive identification of Hemoglobin Pattern


Quick check of pathological samples versus non-pathological samples (color-coded) due to automatic software interpretation

Results examples


FA + Bart's



FS (Sickle Cell Anemia)



FSC (Sickle Cell Disease)



4

Early detection of hemoglobin disorders allows appropriate treatment and follow up, **improving baby's life** while reassuring parents



Ordering information

CAPILLARYS³ DBS

Instruments

- CAPILLARYS 3 DBS1244

Controls/Calibrators

- CAPI 3 NEONAT
Hb AF CONTROL4757
- Hb AFSC CONTROL4792

Reagents


- CAPI 3 NEONAT Hb2506
- CAPI 3 NEONAT HEMOLYSING SOLUTION2563

Consumables

- CAPILLARYS 3 DBS MICROPLATES
(10 units – 96 well barcoded microplates)2578
- CAPILLARYS 3
CAPICLEAN2060
- CAPILLARYS 3
WASH SOLUTION2062
- CAPI 3
REAGENT CUPS (24 x 14)2582



Scan here
for more
information on
CAPILLARYS 3 DBS



For further information, please find your Sebia
local representative to contact on our Website
www.sebia.com/sebia/worldwide-presence/.



Sebia | RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008
Evry cedex | France | Tel.: +33 1 69 89 80 80 | E-mail: sebia@sebia.com | www.sebia.com

©2026 **Sebia** – RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008 Evry cedex | France – Ref: C3DBSBIEN – 04/2023
Non-contractual illustration – Sebia reserves the right to modify at any time and without notice the information contained in this document intended for healthcare
professionals. This document contains information intended for wide distribution and may therefore contain product details or information that is not available or valid in your
country. Carefully read the instructions in the reagent package inserts and instrument manuals.